

Sureseal[®] II

Endoscopic Valve



A large teal geometric shape, consisting of a dark teal triangle and a lighter teal trapezoid, is positioned on the left side of the page, extending from the top left corner towards the bottom right.

THE SEAL THAT KEEPS YOU DRY

The Sureseal II endoscopic valve combines excellent irrigation capability with a leak-tight seal around 0–5F stone baskets, guidewires, laser fibers and other endourological instruments.

Superior Closure

Keeps urologist dry

Works with a variety of scopes and the Applied Forté ureteral access sheaths

Accommodates instruments from 0-5F

Enhanced Handling

Allows passage of instruments into working port without additional manipulation of valve

Enables one-handed manipulation of scope

Latex Free

Protects latex sensitive patients





Sureseal II Endoscopic Valve

Reorder No.	Valve Opening	Pressure Rating	Irrigation Port	Scope Port	Qty/Box
Seal-005	Passive, continuously reliable bidirectional from 0-5F (.065")	20psi or 1034mm Hg	Luer-lock type fitting	Standard luer-lock with ball adapter	10/BOX (includes valve and ball adapter)

Visit [appliedmedical.com/accessories](https://www.appliedmedical.com/accessories) for more information.

The Sureseal II endoscopic valve is used for urological procedures requiring a leak-proof seal.

Please contact your Applied Medical representative for more information on availability. This information is intended for dissemination exclusively to healthcare professionals and is not intended to replace labeling and Instructions for Use (IFU). Please refer to the IFU for the indications, contraindications, warnings, precautions, instructions and other information.

© 2022 Applied Medical Resources Corporation, Rancho Santa Margarita, CA. All rights reserved. Applied Medical, the Applied Medical logo design, and marks designated with the symbol ® are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, the United Kingdom, and/or the European Union.

06/2022 442669-EN-A

