



Epix[®]

Electrosurgical Probes
with Smoke Evacuation

Epix

Electrosurgical Probes with Smoke Evacuation

Probe Tip

Provide precise delivery of monopolar energy




Sliding Sheath

Advances for protection during instrument insertion

Retracts to expose electrosurgical tip


Shafts

Angled and straight variations provide access and improve visibility in small working spaces



Built-In Smoke Evacuation

Evacuates smoke at the push of a button while reducing loss of insufflation



Soft, Flexible Tubing

Drapes easily, reduces coiling and improves control

Recessed Monopolar Electro-surgical Connector

Allows for tangle-free connection to energy generator



Epix Electrosurgical Probe with Smoke Evacuation

Reorder No.	Tip Design	Size	Qty/Box
CW002*	Angled L-Hook	5mm x 42cm	1



Epix Electrosurgical Probe with Smoke Evacuation, Straight

Reorder No.	Tip Design	Size	Qty/Box
CW005**	L-hook	5mm x 42cm	1

* Model CW002 is compatible with Applied Medical GelPOINT® sleeves and Applied Medical Kii® 12 x 60-100mm and 15 x 100-150mm trocars.

** Model CW005 is compatible with Applied Medical GelPOINT sleeves and Applied Medical Kii 5mm trocars (if using non-Kii trocars, a larger port size will be needed).

Visit [appliedmedical.com/epix](https://www.appliedmedical.com/epix) for more information.

Please contact your Applied Medical representative for more information on availability. This information is intended for dissemination exclusively to healthcare professionals and is not intended to replace labeling and Instructions for Use (IFU). Please refer to the IFU for the indications, contraindications, warnings, precautions, instructions and other information.

© 2021 Applied Medical Resources Corporation, Rancho Santa Margarita, CA. All rights reserved. Applied Medical, the Applied Medical logo design, and marks designated with the symbol ® are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, the United Kingdom, and/or the European Union.

12/2021 206425-EN-EUR-C

