



Epix[®]

Laparoscopic
Suction-Irrigation System

Epix

Laparoscopic Suction-Irrigation System

THE ULTIMATE IN COMFORT AND CONTROL

Reusable Probe Attachments



C7202/C7211



C7204

Probe Attachments

Easy to connect and disconnect

Available with a variety of tip offerings and in several lengths and diameters, enabling a wide range of flow rates to be achieved with a single handpiece

Soft, Flexible Tubing

Resistant to coiling and collapsing

Smoke Evacuation at the Source

Easily achieved by simply rotating the red button

Color-Coded Buttons

Seated on soft springs for smooth operation and user control

Clearly differentiated to identify the suction and irrigation functions



Clear, Ergonomic Handpiece

Comfortable to use and transparent, affording an unobstructed view of the working channels



Epix Disposable Tubing Sets

Reorder No.	Description	Diameter	Working Length	Qty/Box
C6001	Probe and valve with dual-bag tubing	5mm	33cm	12
C6003	Valve only with dual-bag tubing*	N/A	N/A	12



Epix Reusable Probes

Reorder No.	Description	Diameter	Working Length	Qty/Box
C7202	Probes with irrigation holes	5mm	33cm	1
C7211	Probes with irrigation holes	5mm	45cm	1
C7204	Probes with irrigation holes	10mm	33cm	1

* For use with reusable probes

Visit [appliedmedical.com/epix](https://www.appliedmedical.com/epix) for more information.

Epix suction-irrigation tubing set and reusable probe attachments are used in procedures to deliver sterile irrigation fluids and to evacuate blood, tissue debris and smoke from the surgical site.

Please contact your Applied Medical representative for more information on availability. This information is intended for dissemination exclusively to healthcare professionals and is not intended to replace labeling and Instructions for Use (IFU). Please refer to the IFU for the indications, contraindications, warnings, precautions, instructions and other information.

© 2022 Applied Medical Resources Corporation, Rancho Santa Margarita, CA. All rights reserved. Applied Medical, the Applied Medical logo design, and marks designated with the symbol ® are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, the United Kingdom, and/or the European Union.

03/2022 310846-EN-EUR-B

