



Forté[®]

Ureteral Access Sheaths

Forté

Ureteral Access Sheaths

Advanced Technology, Clinically Proven

The Forté ureteral access sheaths provide a continuous working channel, simplifying ureteroscopic procedures.

Forté sheaths are clinically proven to protect the urinary system, reduce operating room time, eliminate the need for balloon dilators, and aid in washing out stone fragments.^{1,2,3}



Seamless Transition

The sheath-and-dilator combination offers unprecedented ease of placement.

The dilator ensures easy insertion and reduces trauma.

Gentle Dilation

The tapered dilator gently traverses the ureteral orifice and replaces balloon dilators.

Continual Access

Providing continual access through a sheath reduces procedural time.

Using a sheath protects the scope and instruments, reducing the need for costly repairs.

Less Trauma

The risk of trauma decreases when the ureter is protected by a sheath during ureteroscopy.

The hydrophilic coating facilitates advancement.



Forté Ureteral Access Sheaths

Reorder No.	Dilator	Sheath	Length	Qty/Box
<i>AxP</i>				
B7276	10F	12F-16F	28cm	6
B7216	10F	12F-16F	35cm	6
B7296	10F	12F-16F	45cm	6
B7286	10F	12F-16F	55cm	6



Reorder No.	Dilator	Sheath	Length	Qty/Box
<i>AxP</i>				
B7066	12F	14F-18F	20cm	6
B7016	12F	14F-18F	35cm	6
B7096	12F	14F-18F	45cm	6
B7086	12F	14F-18F	55cm	6



Reorder No.	Dilator	Sheath	Length	Qty/Box
<i>AxP</i>				
B7166	14F	16F-18F	20cm	6
B7176	14F	16F-18F	28cm	6
B7116	14F	16F-18F	35cm	6

1. Auge BK, Pietrow PK, Lallas CD, Raj GV, Santa-Cruz RW, Preminger GM. Ureteral access sheath provides protection against elevated renal pressures during routine flexible ureteroscopic stone manipulation. *Journal of Endourology*. 2004;18(1):33-36. doi:10.1089/089277904322836631.
2. Kourambas J, Byrne RR, Preminger GM. Does a ureteral access sheath facilitate ureteroscopy?. *Journal of Urology*. 2001;165(3):789-793.
3. Ajayi I, Lillas L, Alan J, Montgomery B, Naerger H, Palfrey E. The use of ureteral access sheaths and flexible ureteroscopy; a marriage in heaven?[D3-P7] *Journal of Endourology*. 2001;15(suppl):S1.

Visit [appliedmedical.com/forte](https://www.appliedmedical.com/forte) for more information.

The Forté ureteral access sheaths provide a working channel in ureteroscopic procedures.

Please contact your Applied Medical representative for more information on availability. This information is intended for dissemination exclusively to healthcare professionals and is not intended to replace labeling and Instructions for Use (IFU). Please refer to the IFU for the indications, contraindications, warnings, precautions, instructions and other information.

© 2022 Applied Medical Resources Corporation, Rancho Santa Margarita, CA. All rights reserved. Applied Medical, the Applied Medical logo design, and marks designated with the symbol ® are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, the United Kingdom, and/or the European Union.