



# GelPort<sup>®</sup>

Laparoscopic System

# GelPort

## Laparoscopic System

---

### GelSeal® Technology

Enables unlimited hand and instrument exchanges while maintaining pneumoperitoneum.

Accommodates trocar insertion for enhanced versatility and patient cosmesis.

Features a small profile for limited abdominal space.



### Enabling Design

Offers unparalleled access for rapid dissection and mobilization of tissue.

Facilitates a Pfannensteil approach to optimize advanced procedures.

Facilitates simple and rapid setup.



### Alexis® Wound Protector-Retractor

Adjusts to retract varying abdominal walls and different incision sizes.

Offers 360° of atraumatic retraction and protection<sup>1-2</sup> for enhanced exposure and access.

Simplifies extracorporeal resection and specimen retrieval.

Allows clear visualization of wound margins.

1. Reid K., et al. Barrier Wound Protection Decreases Surgical Site Infection in Open Elective Colorectal Surgery: A Randomized Clinical Trial. *Dis Colon Rectum*. 2010 Oct; 53(10): 1374-1380.

2. Horiuchi T., et al. Randomized Controlled Investigation of the Anti-Infective Properties of the Alexis Retractor/Protector of Incision Sites. *J Trauma*. 2007 Jan; 62(1): 212-215.



### **Union of Exceptional Technologies**

Allows unrestricted access for surgical speed and efficiency.

Simplifies specimen retrieval.

Provides 360° of atraumatic retraction and protection.<sup>1-2</sup>

Offers trocar access for unparalleled versatility.



## GelPort Laparoscopic System


Reorder No.	Qty/Box
-------------	---------

C8XX2CE	1
---------	---

### *Components*

(1) GelSeal Cap
-----------------

(1) Alexis wound protector-retractor
--------------------------------------



Visit [www.appliedmedical.com/gelport](http://www.appliedmedical.com/gelport)  
for more information

Devices listed may not be approved in all markets. Please contact your Field  
Implementation team member for more information on availability.

© 2019 Applied Medical Resources Corporation. All rights reserved. Applied Medical, the Applied  
Medical logo design and marks designated with a ® are trademarks of Applied Medical Resources  
Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South  
Korea, the United States and/or the European Union  
227644-EN-EUR-A