



Applied Scope Warmer

A large teal geometric shape, consisting of a dark teal triangle and a lighter teal trapezoid, is positioned on the left side of the page, extending from the top-left corner towards the bottom-right.

A SIMPLE SOLUTION FOR BETTER VISIBILITY

Simple to use
Minimises lens fogging

The Answer to Laparoscopic Fogging

Eliminates the need for chemical gels or wipes, retains heat throughout lengthy procedures and minimises the need for repeated scope cleansing and defogging.

Scope Warmer Seal

Retains heat throughout lengthy procedures

Scope Warmer Bottle

Minimises fogging by pre-warming the laparoscope to body temperature

Scope Warmer Base

Stabilises the system when a scope is inserted





Applied Scope Warmer

Reorder No.	Description	Qty/Box
C3001	Scope Warmer Bottle* Reusable Volume Capacity = 500mL	1 EACH
C3002	Scope Warmer Base* Reusable	1 EACH
C3101	Scope Warmer Seals Disposable, Single Use Individually packaged in sterile pouches	20

*Recommended Sterilization Method: Autoclave or STERRAD** 1100S system.
Please see Instructions for Use (IFU) for complete information.
Scope is not included.

Visit [appliedmedical.com/laparoscopic](https://www.appliedmedical.com/laparoscopic) for more information.

Please contact your Applied Medical representative for more information on availability. This information is intended for dissemination exclusively to healthcare professionals and is not intended to replace labeling and Instructions for Use (IFU). Please refer to the IFU for the indications, contraindications, warnings, precautions, instructions and other information.

© 2021 Applied Medical Resources Corporation. All rights reserved. Applied Medical, the Applied Medical logo design, and marks designated with the symbol © are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, the United Kingdom, and/or the European Union.

12/2021 414709-EN-EUR-A

