



Insufflation Needle

A large teal graphic element consisting of a diagonal line that starts from the top left and extends towards the bottom right, creating a triangular shape on the left side of the page. The teal color is a medium-dark shade.

ENABLES HIGH FLOW INSUFFLATION

13 gauge percutaneous needle

Tactile and visual indicators

Ergonomic housing design

Insufflation Needle



Spring Loaded Blunt Stylet with Red Indicator Button

Provides tactile and visual feedback when stylet is retracted

13 Gauge High Flow Needle

Offers faster insufflation than traditional needles, reducing OR time



Ergonomic Housing Design

Ensures control during needle insertion





Insufflation Needle

Recorder No.	Size	Quantity/Box
C2201	13 Gauge, 120mm	20
C2202	13 Gauge, 150mm	20

Visit [appliedmedical.com/Laparoscopic](https://www.appliedmedical.com/Laparoscopic) for more information.

Please contact your Applied Medical representative for more information on availability. This information is intended for dissemination exclusively to healthcare professionals and is not intended to replace labeling and Instructions for Use (IFU). Please refer to the IFU for the indications, contraindications, warnings, precautions, instructions and other information.

© 2021 Applied Medical Resources Corporation. All rights reserved. Applied Medical, the Applied Medical logo design, and marks designated with the symbol © are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, the United Kingdom, and/or the European Union.

11/2021 411930-EN-EUR-A

